

# Saving a life deserves the best protection



## Defib Store 1000 Series



The Defib Store 1000 is a robust AED cabinet at an affordable price. Made from high grade materials by a leading cabinet specialist.

- Outdoor defibrillator cabinet
- Available in mild steel or stainless steel
- Completely water and dust tight IP66
- Vandal proof features
- Marine grade stainless steel keypad lock
- Thermostatically controlled heater
- Motion sensor LED light
- 10 Year guarantee on Stainless Steel model
- Made in the UK



[www.defibstore.co.uk](http://www.defibstore.co.uk)

# Defib Store recognise the importance of a high-quality, secure cabinet to store a public access defibrillator

Part of the CE-TEK group, a specialist enclosure manufacturer for over 30 years, Defib Store has utilised CE-TEK's expertise to develop a range of outdoor, robust defibrillator cabinets for all environments. Our industry leading range of weatherproof cabinets ensure defibrillators can be secure and always accessible

The popular Defib Store 1000 series is now available in mild or stainless steel

Defib Store 1000 features include:


Water and dust tight. IP66 rated to withstand even the harshest of weather conditions

Thermostatically controlled heater to ensure the AED is always within it's working temperature range


Marine grade stainless steel keypad lock; tested with salt spray for 1,000 hours, ensures the lock remains fully functional in all weather conditions

Battery powered motion sensor LED light to illuminate cabinet on opening. Easily detachable for use as a torch in a rescue


Secure viewing window ensures AED status indicator is visible externally at all times



Defib Store 1000 lle



Defib Store 1000 ule



Defib Store 1000 SS lle

Defib Store 1000 SS ule

Internal AED hook accommodates any defibrillator (without carry case)

Available unlocked with a plastic winged handle and optional alarm\*

Powder coated paint in high visibility yellow or green as standard. Also available in any paint colour\*

All Defib Store cabinets include printed graphics on high grade laminated vinyl with the option to customise to include personalised logos/text as required\*

Cabinets are delivered fully assembled with UK plug, fixing kit, installation and maintenance guide

2 year guarantee on mild steel cabinet, 10 year guarantee on stainless steel cabinet


\* Small additional charge

Cabinet External Dimensions: 400mm (H) x 400mm (W) x 200mm (D)

Product code	Material	Colour	Door Locking	Heating & Lighting	IP Rating	Viewing Window	Location
DFS1K442LLE	Mild Steel	Yellow/Green	SS Keypad lock	Yes	IP66	Yes	Indoor/Outdoor
DFS1K442ULE	Mild Steel	Yellow/Green	Plastic Wing Handle	Yes	IP66	Yes	Indoor/Outdoor
DFS1K442LLE-SS	Stainless Steel	Yellow/Green	SS Keypad lock	Yes	IP66	Yes	Indoor/Outdoor
DFS1K442ULE-SS	Stainless Steel	Yellow/Green	Plastic Wing Handle	Yes	IP66	Yes	Indoor/Outdoor

All Defib Store 1000 models available without electrics (heating & lighting)

Please note Defibrillators installed within 5 miles from the coast should be installed in our stainless steel cabinets



[www.defibstore.co.uk](http://www.defibstore.co.uk)

Defib Store Ltd, Units 1-4 Tideswell Business Park, Tideswell, Derbyshire SK17 8NY

Tel: +44 (0) 1298 872186

Fax: +44 (0) 1298 871155

E-mail: [sales@defibstore.co.uk](mailto:sales@defibstore.co.uk)